



SOOO ADHERENT IS MUCOPREN® SOFT

Holds to its promises

Mucopren® Soft is a soft long-term relining material for removable dentures based on A-silicone for direct and indirect application. The material is also very useful for many other indications, like relieving the operation site after placing implants, for compensating areas of mucous membrane with varying degrees of resiliency and also for relieving sharp bony ridges.



User advantages


- Resistant to contamination*
- Can be precisely trimmed
- Excellent bond strength to the denture
- For direct and indirect application
- Simple handling

Product properties

- Excellent adhesion to acrylic dentures through patented Mucopren® Adhesive
- Smooth, hydrophobic silicone surface*
- Silicone-based sealing
- High tear-resistance

Benefits for the dentist

- Optimum requirements for hygienic denture conditions
- Long period in situ and comfort
- Modern material characteristics ensure easy, reliable handling
- Very efficient, value creation remains in the dental practice



For denture relining
without compromises!

Benefits of using the new adhesive

ONLY A SINGLE
APPLICATION
NECESSARY



CHAIRSIDE

3 MIN

The processing time is reduced to only **3 minutes** following mouth removal. No water bath necessary.



LABSIDE

3 MIN

The processing time is reduced to only **3 minutes** in the water bath or pressure pot.



* Concluded from in vitro studies which showed that the smooth silicone surface prevented Candida albicans growth.